## **SPECIFICATIONS**

- > Pre-assembled Electromyography (EMG) sensor with **SnapBIT-DUO**
- > Fixed electrode distance
- > Single cable connection

#### **FEATURES**

- > Allows basic EMG data acquisition
- > Easy to wear
- > Plug & play design
- > Raw data acquisition

## **APPLICATIONS**

- > Human-Computer Interaction
- > Robotics & Cybernetics
- > Physiology studies
- > Psychophysiology
- > Biomechanics
- > Biofeedback
- > Muscle reflex studies
- > Nerve conduction measurement
- > Biomedical devices prototyping

### GENERAL DESCRIPTION

The BITalino assembled EMG sensor is designed for everyone who wants to measure muscle activity by evaluating electromyography (EMG) signals. This bundle is completely assembled with our 3D Printed Casing for BITalino (r)evolution Plugged making it more convenient to use, wearable, sharable & transportable. The Assembled Electromyography (EMG) SnapBIT-DUO with Sensors allow repeatedly accurate & fast measurements, once the user can benefit from the prefixed electrode distance, 1-lead cable for each sensor (instead of 2-lead or 3-lead cable). Only one reference cable is needed when working with one or several assembled sensors.



Fig. 1. Assembled EMG sensor - top view.



Fig. 2. Assembled EMG sensor – top (left) and bottom (right) view from the inside.




Fig. 3. Typical raw EMG data (acquired with BITalino (r)evolution) on the muscle biceps brachii with flexion of elbow and forearm in supination.



Fig. 4. Example electrode placement of the assembled version along the muscle fibers and the reference electrode on a bone.



PLUX - Wireless Biosignals, S.A. Av. 5 de Outubro, n. 70 - 2 1050-059 Lisbon, Portugal bitalino@plux.info http://bitalino.com/

**REV A** 

© 2020 PLUX 😉 🕦



This information is provided "as is," and we make no express or implied warranties whatsoever with respect to functionality, operability, use, fitness for a particular purpose, or infringement of rights. We expressly disclaim any liability whatsoever for any direct, indirect, consequential, incidental or special damages, including, without limitation, lost revenues, lost profits, losses resulting from business interruption or loss of data, regardless of the form of action or legal theory under which the liability may be asserted, even if advised of the possibility of such damages.





# **Electromyography (EMG) Assembled Sensor Data Sheet**

## TRANSFER FUNCTION

**EMG** datasheet

# PHYSICAL CHARACTERISTICS

**EMG** datasheet

PACKAGING Weight: 15 g

## ORDERING GUIDE

Part #	Description
BUNDLE-EMG-	Pre-assembled sensor for muscle activity measurement.
UCE6	·

