

Electrocardiography (ECG) Sensor Datasheet

ECG 10082020

SPECIFICATIONS

- > **Gain:** 1019
- > **Range:** $\pm 1.47\text{mV}$ (with $V_{CC}=3\text{V}$)
- > **Bandwidth:** 25-100Hz
- > **Input Impedance:** $>100\text{G}\Omega$
- > **CMRR:** 100dB
- > **Cable Length:** $100\pm 0.5\text{cm}$ (customizable)
- > **Connector Type:** UC-E6 (Male)

FEATURES

- > Bipolar differential measurement
- > Pre-conditioned analog output
- > High signal-to-noise ratio
- > Medical-grade raw data output
- > Ready-to-use & miniaturized form-factor

APPLICATIONS

- > Life sciences studies
- > Heart rate & heart rate variability
- > Human-Computer Interaction
- > Biometrics
- > Affective computing
- > Physiology studies
- > Psychophysiology
- > Biofeedback
- > Biomedical devices prototyping

GENERAL DESCRIPTION

Our low-noise ECG local differential triode configuration enables fast application and unobtrusive single-lead ECG data acquisition (although custom electrode cable configurations are available). The state-of-the-art design of the analog frontend on this sensor is specifically targeted at analyzing minutiae in the data and provides medical-grade raw sensor data.

This sensor can be used to extract heart rate data and other ECG features, enabling its application in research fields such as biomedical, biofeedback, psychophysiology, and sports, among many others.




Fig. 1. Standard biosignalsplux ECG sensor with short electrode cables (4cm + 6cm + 4cm).




Fig. 2. Typical raw ECG data (acquired with biosignals).




Fig. 3. Example sensor placement using an ECG triode (equivalent to a standard medical-grade V6 lead).

biosignalsplux
wearable body sensing platform

PLUX – Wireless Biosignals, S.A.
Av. 5 de Outubro, n. 70 – 2.
1050-059 Lisbon, Portugal
plux@plux.info
<http://biosignalsplux.com/>

REV C

© 2020 PLUX

This information is provided "as is," and we make no express or implied warranties whatsoever with respect to functionality, operability, use, fitness for a particular purpose, or infringement of rights. We expressly disclaim any liability whatsoever for any direct, indirect, consequential, incidental or special damages, including, without limitation, lost revenues, lost profits, losses resulting from business interruption or loss of data, regardless of the form of action or legal theory under which the liability may be asserted, even if advised of the possibility of such damages.

Electrocardiography (ECG) Sensor Datasheet

ECG 10082020

TRANSFER FUNCTION

[-1.47mV, +1.47mV]

$$ECG(V) = \frac{\left(\frac{ADC}{2^n} - \frac{1}{2}\right) \times VCC}{G_{ECG}}$$

$$ECG(mV) = ECG(V) \times 1000$$

$VCC = 3V$ (operating voltage)

$G_{ECG} = 1019$ (sensor gain)

$ECG(V)$ – ECG value in Volt (V)

$ECG(mV)$ – ECG value in millivolt (mV)

ADC – Value sampled from the channel

n – Number of bits of the channel¹

ELECTRODE CONNECTIONS

Sleeve Color	Red	Black	White
Electrode Cable	+	-	Reference

PHYSICAL CHARACTERISTICS

> $W_1 \times L_1 \times H_1$: 1.5 x 2.1 x 0.4cm

> $W_2 \times L_2 \times H_2$: 1.4 x 1.7 x 0.5cm

> A_1 : 105±0.5cm


> A_2 : 4±0.5cm

> A_3 : 6±0.5cm

> D : 0.4cm

> **Weight:** 9g

> **Available Sleeve Colors:** White, Black, Blue, Green, Red, Yellow, Gray, and Brown



Note: This sensor is also available with the electrodes cables (A_2 & A_3) in 30cm lengths.

¹ The number of bits for each channel depends on the resolution of the Analog-to-Digital Converter (ADC); in biosignalsplux the default is 16-bit resolution ($n = 16$), although 12-bit ($n = 12$) and 8-bit ($n = 8$) may also be found on older devices.

Electrocardiography (ECG) Sensor Datasheet

ECG 10082020

ORDERING GUIDE

Reference	Package Description
SENSPRO-ECG1	Electrocardiography (ECG) sensor with standard physical characteristics and a random cable sleeve color

Electrocardiography ECG (3 x 30 cm)

Sensor Data Sheet

ECG 26012018

SPECIFICATIONS

- >Gain: 1000
- >Range: $\pm 1.5\text{mV}$ (with VCC = 3V)
- >Bandwidth: 0.5-100Hz
- >Consumption: $\sim 1\text{mA}$
- >Input Impedance: $>100\text{G}\Omega$
- >CMRR: 100dB

FEATURES

- >Bipolar differential measurement
- >Pre-conditioned analog output
- >High signal-to-noise ratio
- >Shielded miniaturized cables
- >Medical-grade raw data output
- >Ready-to-use form factor

APPLICATIONS

- > Life sciences studies
- > Heart rate & heart rate variability
- > Human-Computer Interaction
- > Biometrics
- > Affective computing
- > Physiology studies
- > Psychophysiology
- > Biofeedback
- > Biomedical devices prototyping

GENERAL DESCRIPTION

Electrocardiography (ECG) records electrical activity of the heart over time. Variations in the duration, amplitude, and morphology of the ECG waves are used for diagnosing abnormal cardiac rhythms and conduction patterns. Our low-noise ECG local differential triode configuration enables fast application and unobtrusive data acquisition (although custom electrode cable configurations are available). The state-of-the-art design of the analog frontend on this sensor is specifically targeted at analyzing minutiae in the data. Together with the Heart Rate Variability (HRV) plugin on our OpenSignals software, one can easily record and extract meaningful information from the collected data. Examples:

<http://bit.ly/1ddQnsv>
<http://bit.ly/1JEW2lk>




Fig. 1. Triode electrode configuration for fast, minimally intrusive setup on your subjects.




Fig. 2. Typical raw ECG data (acquired with biosignals).




Fig. 3. Example sensor placement (equivalent to a standard medical-grade lead I).

biosignalsplux
wearable body sensing platform

PLUX – Wireless Biosignals, S.A.
Av. 5 de Outubro, n. 70 – 8.
1050-059 Lisbon, Portugal
plux@plux.info
<http://biosignalsplux.com/>

© 2018 PLUX

This information is provided "as is," and we make no express or implied warranties whatsoever with respect to functionality, operability, use, fitness for a particular purpose, or infringement of rights. We expressly disclaim any liability whatsoever for any direct, indirect, consequential, incidental or special damages, including, without limitation, lost revenues, lost profits, losses resulting from business interruption or loss of data, regardless of the form of action or legal theory under which the liability may be asserted, even if advised of the possibility of such damages.

Electrocardiography ECG (3 x 30 cm) Sensor Data Sheet

TRANSFER FUNCTION

[-1.5mV, 1.5mV]

$$ECG(V) = \frac{\left(\frac{ADC}{2^n} - \frac{1}{2}\right) \times VCC}{G_{ECG}}$$

$$ECG(mV) = ECG(V) \times 1000$$

$VCC = 3V$ (operating voltage)

$G_{ECG} = 1000$ (sensor gain)

$ECG(V)$ – ECG value in Volt (V)

$ECG(mV)$ – ECG value in millivolt (mV)

ADC – Value sampled from the channel

n – Number of bits of the channel¹

PHYSICAL CHARACTERISTICS

> **W1 x L1 x H1**: 1.0x1.8x0.4cm

> **W2 x L2 x H2**: 1.5x2.3x0.4cm

> **A1**: 100.0±0.5cm


> **A2**: 30.0±0.5cm

> **A3**: 30.0±0.5cm

> **A4**: 30.0±0.5cm

> **D**: 0.4cm

> **S**: White, Black, Blue, Green, Red, Yellow, Gray, or Brown



¹The number of bits for each channel depends on the resolution of the Analog-to-Digital Converter (ADC); in biosignalsplux the default is 16-bit resolution ($n = 16$), although 12-bit ($n = 12$) and 8-bit ($n = 8$) may also be found.

Electrocardiography ECG (3 x 30 cm) Sensor Data Sheet

ORDERING GUIDE

Reference	Package Description
SENSPRO-ECG1	Electrocardiography (ECG) sensor with standard physical characteristics and a random cable sleeve color.
SENSPROECG1-100-30-30-S	Electrocardiography (ECG) sensor built with electrode cables of 30cm (A2 & A3 in illustration above) and custom sleeve color.
SENSPROECG1-A1-A2-A3-S	Electrocardiography (ECG) sensor built with custom lengths for A1, A2, and A3 (all in cm), and custom sleeve color S; for standard physical characteristics in A1, A2, A3, or S use 0. Examples: > ECG1-200-0-0-0: Otherwise all-standard ECG sensor except for a 200cm cable A1 > ECG1-0-0-0-Yellow: Otherwise all-standard ECG sensor except for a yellow cable sleeve > ECG1-50-10-10-Red: Fully custom ECG sensor with a 50cm cable A1, 10cm electrode cables A2 & A3, and a red cable sleeve