

Acoustic Sensor Data Sheet

ACOUSTIC 21032018

SPECIFICATIONS

- > **Type:** High-performance, low-power, silicon microphone
- > **Range:** 0-3V (0-1 normalized)
- > **Gain:** 410
- > **Bandwidth:** 29Hz – 1591Hz

FEATURES

- > Audio power measurement
- > Silicon microphone
- > Pre-conditioned analog output
- > Small form-factor
- > Raw data output

APPLICATIONS

- > Biomedical research
- > Auditive gastric activity monitoring
- > Audio power measurement
- > Sound (impulse) detection

GENERAL DESCRIPTION

The biosignalsplux Acoustic Sensor is designed for audio power and audio impulse measurements, which can be used to monitor ambient noises, physiological sounds, and/or sound impulses. This sensor can be used in combination with the biosignalsplux acquisition unit and other biosignalsplux sensors (e.g. ECG), to track the changes of physiological signals as results of induced auditive stimuli or changing activity of physiological processes (e.g. abdominal sounds caused by gastrointestinal activity).




Fig. 1. Typical raw audio data of a series of sound impulses (acquired with biosignalsplux and the acoustic sensor).

biosignalsplux
wearable body sensing platForm

PLUX – Wireless Biosignals, S.A.
Av. 5 de Outubro, n. 70 – 8.
1050-059 Lisbon, Portugal
plux@plux.info
<http://biosignalsplux.com/>

REV A

© 2018 PLUX

This information is provided "as is," and we make no express or implied warranties whatsoever with respect to functionality, operability, use, fitness for a particular purpose, or infringement of rights. We expressly disclaim any liability whatsoever for any direct, indirect, consequential, incidental or special damages, including, without limitation, lost revenues, lost profits, losses resulting from business interruption or loss of data, regardless of the form of action or legal theory under which the liability may be asserted, even if advised of the possibility of such damages.

SENSOR (SEN)

Sensor Data Sheet

SEN XXXX2018

TRANSFER FUNCTION

[0V, 3V]

$$ACOUSTIC(V) = \frac{\left(\frac{ADC}{2^n - 1} - \frac{1}{2}\right) \times VCC}{G_{Acoustic}}$$

Normalized: [0, 1]

$$ACOUSTIC(-) = \frac{ADC}{(2^n - 1)}$$

- ACOUSTIC(V)* – Sensor output in V
- ACOUSTIC(-)* – Normalized sensor output
- ADC* – Value sampled from the channel
- n* – Number of bits of the channel¹
- VCC* – 3V (operating voltage)

PHYSICAL CHARACTERISTICS

- > **Size (W x L x H):** 1.0x1.8x0.4cm
- > **Cable Length (A):** 105.0±0.5cm
- > **Weight:** 9g±1g
- > **Sleeve:** White, Black, Blue, Green, Red, Yellow, Gray, or Brow




Fig. 2. Physical characteristics of the sound sensor.

ORDERING GUIDE

Reference	Package Description
SENSPRO-ACOUSTIC-1	Acoustic sensor with standard physical characteristics and random sleeve color.
SENSACOUSTIC-1-A-S	Acoustic sensor with custom cable length A and/or custom sleeve color S; for standard physical characteristics for A or S use 0.

¹ The number of bits for each channel depends on the resolution of the Analog-to-Digital Converter (ADC); in biosignalsplux the default is 16-bit resolution ($n = 16$), although 12-bit ($n = 12$) and 8-bit ($n = 8$) may also be found.