

Temperature (TMP) Sensor Data Sheet

TMP 200616

SPECIFICATIONS

- > **Range:** [-40°C; +125°C]
- > **Accuracy:** $\pm 2^\circ\text{C}$
- > **Linearity:** $\pm 0.5^\circ\text{C}$
- > **Consumption:** $\sim 0.05\text{mA}$

FEATURES

- > Calibrated directly in $^\circ\text{C}$
- > Small form factor
- > Raw data output
- > Easy-to-use

APPLICATIONS

- > Body temperature measurement
- > Environmental analysis
- > Psychophysiology
- > Biofeedback
- > Biomedical devices prototyping

GENERAL DESCRIPTION

Our calibrated Temperature (TMP) sensor has been especially selected for seamless operation with BITalino (r)evolution, targeting the measurement of body or environmental temperature. Its small form factor enables easy application on any surface. Typical applications include peripheral temperature data acquisition.




Fig. 1. Pin-out and physical dimensions.

bitalino

PLUX – Wireless Biosignals, S.A.
Av. 5 de Outubro, n. 70 – 8.
1050-059 Lisbon, Portugal
bitalino@plux.info
<http://bitalino.com/>

REV A

© 2016 PLUX 

This information is provided "as is," and we make no express or implied warranties whatsoever with respect to functionality, operability, use, fitness for a particular purpose, or infringement of rights. We expressly disclaim any liability whatsoever for any direct, indirect, consequential, incidental or special damages, including, without limitation, lost revenues, lost profits, losses resulting from business interruption or loss of data, regardless of the form of action or legal theory under which the liability may be asserted, even if advised of the possibility of such damages.



BEWARE: DIRECT OR INDIRECT COUPLING TO THE MAINS MAY RESULT IN SHOCKING HAZARD



Temperature (TMP) Sensor Data Sheet

TRANSFER FUNCTION

[-40°C, +125°C]

$$TMP(^{\circ}C) = \left(\frac{ADC}{2^n} \cdot VCC - 0.5 \right) \cdot 100$$

$$TMP(^{\circ}F) = TMP(^{\circ}C) \cdot \frac{9}{5} + 32$$

$VCC = 3.3V$ (operating voltage)

$TMP(^{\circ}C)$ – TMP value in degrees Celsius ($^{\circ}C$)

$TMP(^{\circ}F)$ – TMP value in degrees Fahrenheit ($^{\circ}F$)

ADC – Value sampled from the channel

n – Number of bits of the channel¹

ORDERING GUIDE

Part #	Description
SENS-TMP-NC	Temperature (TMP) sensor without connectors
SENS-TMP-UCE6	Temperature (TMP) sensor with UC-E6 socket for seamless plug & play connection to a BITalino (r)evolution Plugged or Core
SENS-TMP-SHER	Temperature (TMP) sensor with a Molex Sherlock 4-pin socket for easy power and signal cable connection or pin breakout using PCB wires

¹ The number of bits for each channel depends on the resolution of the Analog-to-Digital Converter (ADC); in BITalino the first four channels are sampled using 10-bit resolution ($n = 10$), while the last two may be sampled using 6-bit ($n = 6$).